

## Hose Assemblies For Medical Gases



Diss wrench style  
nut & nipple



Ohmeda style female  
quick connect



Allied style male  
Adaptor



Ohmeda style male  
quick connect



DISS quick  
connect



Hand tight DISS  
nut & nipple




Puritan style  
male adaptor

**Size of Hose** (all hose assemblies listed are conductive except Evacuation is non-conductive)

4 = 1/4" hose

90 = 90 degree angle

5 = 5/16" hose



### Gas Service

04 = Nitrous Oxide  
08 = Carbon Dioxide  
12 = Nitrogen  
16 = Medical Air  
22 = Vacuum  
24 = Oxygen  
32 = Evacuation

### Please use the following codes for both inlet and outlet:

MM = Ohmeda style male adaptor  
AM = Allied Chemetron style male adaptor  
PM = Puritan style male adaptor  
EM = Oxequip "07" style male adaptor  
RM = Oxequip "Medstar" style male adaptor  
SM = Schrader swivel style male adaptor  
MF = Ohmeda style female coupler  
AF = Allied Chemetron style female coupler  
PF = Puritan style female coupler  
SF = Schrader style female coupler  
DV = DISS male demand valve (one check)  
BA = DISS male body adaptor  
WT = DISS wrench tight female nut & nipple  
HT = DISS hand tight nut & nipple  
ZN = DISS quick connect hand tight nut & nipple  
2M = 1/8" male pipe thread  
4M = 1/4" male pipe thread  
2F = 1/8" female pipe thread  
4F = 1/4" female pipe thread  
XX = blank end

### Example:

Part # H424096WTMMR  
This is a 1/4" ID Green OXY hose with a DISS 1240 wrench style nut and nipple on one end and an Ohmeda style male on the other end and a double 'key chain' type hose retractor attached to the hose